

Area	Sub-Area	Device	Level of Disinfection	Freq	Who	Agent
Professional	Exam Room	Spuds, Tweezers, Alger Brush	High	P	O	B10
		Tonometer Pachimeter probe	High	P	O	B10
		Gonioscopy Lenses	High	P	O	B10
		Contact Fundus Lenses	High	P	O	B10
		Forehead/chinrest OCT, camera, NCT, AR	Low	P	O/S	A
		Occluders, eye patch	Low	P	O	B10
		Diagnostic Equipment (perimeter, OCT)	Low	P	O/S	B50
		Sinks	Low	2	O/S	B50
		Exam Chair/Unit	Low	P	O	B50
		R/G glasses, polar glasses	Low	P	O	B50
		Trial frame/lenses	Low	P	O	B50
		Hand held instruments	Low	2	O	B50
		78/90/20D	Low	P	O	B50
	Lab/Disp Area					
		Frame warmer	Low	2	S	B50
		Frames on display	Low	P	S	S
		Frame Displays	Low	2	S	B50
		Lab hand tools	Low	2	S	B50
Administrative	Desk Counters		Low	2	S	B50
	Computer Keyboards, Mouse, Telephone		Low	2	S	B50
	Visa Device		Low	P	S	L
	Staplers, Tape Dispenser		Low	2	S	B50
	Pens		Low	2	S	B50
	Fax Machine		Low	2	S	B50
General Office	Waiting Area		Low	2	S	B50
	Door Handles		Low	2	S	B50
	Washrooms		Low	2	S	F
	Light Switches		Low	2	S	B50

Frequency Legend:
P: After direct patient contact
2: Twice a day

Disinfection Agent:
B10: 1:10 Bleach
B50: 1:50 Bleach
F: Fantastik
Q: Quat solution
A: Alcohol
S: Dish Soap
L: Lysol wipe

Who Legend:
O: Optometrist
S: Staff